

Lÿnfinity

Spray Freeze Drying

Presenter: Dennis Posheluk, Technical Product Manager, IMA Life



November 2020

IMA  **LIFE**

Aseptic Processing & Freeze Drying Solutions


- The design of freeze dryers has not changed much in the last 30 years
 - Significant changes around steam sterilization of the equipment and clean-in-place (CIP/SIP)
 - Use of non-CFC based refrigerants (for global warming potential); reduced thermal mass from stainless steel and vessel construction
 - Auto-loading system and operator/process risk mitigation with the use of RABS and followed by barrier-isolator technology
 - Improved use of PAT for real-time process analysis
- All of the above developments are aimed at improving/maintaining the product quality/operator safety

Traditional Freeze drying batch process = time consuming and energy intensive

- Drying rate is limited by poor heat transfer coefficient leading to long cycles.
- While the use of new process analytical technology has aided in-process understanding, most processes are run conservatively, making it inherently inefficient.

Thus, there is a need to re-think the heat and mass transfer for making the process more efficient.

SCHEMATIC: FREEZE-DRYING



Drying is limited by the mass transfer resistance offered by the drying cake

Continuous Spray Freeze Drying process = high throughput and high cycle efficiency

The current development effort focuses on developing a robust yet gentle continuous aseptic process for spray freeze-drying.

Spray Freeze Drying (SFD) can have a transcending impact on process efficiency related to Biologic drug substance storage, inhalation systems, antibiotics for bulk storage among others with direct powder dosing capabilities directly into product containment systems.



Vials



Dual chamber cartridge



Inhalation device



Dual chamber syringe

Trend in container systems

- Antibiotics and other Bulk dried products require milling and powder filling today.
- Better self administration through the use of pre-filled syringe/cartridge system with improved process efficiency.
- Most Biologic drug substance (DS) are often stored and transported in cryogen vessels (>100 L) before being converted to drugs, making it extremely expensive.
- Vaccines require high volume, low downtime/cost production
- Opportunity to reduce recon time in hard to reconstitute products.



- It is no longer necessary to ship and store freeze-dried vaccines (measles, yellow fever and BCG) at -20°C. Instead, they may be refrigerated at +2° to +8°C.



**DEPARTMENT OF VACCINES
AND BIOLOGICALS**



World Health Organization
Geneva
2000




Freezing:

Spray of liquid formulation into a cryogenically cooled freezing chamber.

Drying:

The resultant frozen particles are freeze-dried in the drying chamber.

The drying chamber contents are both transported and heated in a vacuum chamber to promote rapid sublimation and prevent agglomeration.



LYNFINITY


A continuous spray freeze drying technology that produces a dry powder in bulk form aseptically.

Technology

- Spray freezing in a freezing chamber.
- Continuous freeze drying of frozen particles in vibratory agitated drying chamber.

Resulting dry powder is ready for dosing into vial, syringe or inhaler.

Possibility to use CIP/SIP powder filler.




Spraying Step


- The product spray is initiated when a steady laminar product feed is broken into uniform droplets at the top of the freezing chamber.

Freezing Step

- Introduction of spray product into LN₂ cooled freezing tower.
- Product particles are quickly frozen and accumulated at bottom of container.
- Particle size is controlled for uniformly frozen particles, allowing rapid drying in subsequent steps.



Freezing Chamber




High speed camera view of frozen particles



Drying Step

- Frozen particles are transferred into a temperature and vacuum controlled drying chamber.
- Surface heating is provided by individually temperature controlled shelves.
- Transport of the particles is done by vibrating shelves.
- Dual condenser design to allow continuous operation with conventional vacuum design.





- Fills sterile powder in low cost Grade C/D background
- CIP-SIP compatible

- This laboratory module for spray freezing features a controlled spray of liquid product into a cryogenically cooled freezing column
- The design of the laboratory freezing module is identical to that used in the production version to allow for easy model building and cycle development
- The frozen pellets are manually transferred into a separate drying module such as a traditional laboratory freeze dryer



The LYnfinity LAB offers access to the emerging technology of spray freeze drying:

- Ease of use
- Low barrier to entry prior to moving product towards continuous freeze drying
- Generate uniform frozen droplets of mean volumetric diameter 600 microns
- Use of traditional freeze dryer as drying module
- Process development with limited formulation (down to 200 ml)
- Technology evaluation
- Model building
- Dried product characterization
- Stability testing
- Small footprint for setup in the lab
- Designed for handle and install in laboratory with standard double doors





(1)



(2)



(3)

- (1) User loads liquid product into liquid product reservoir
- (2) User initiates LYnfinity LAB spray freeze cycle through HMI panel, with cycle status displayed through the HMI and LED strip
- (3) Upon cycle completion, user removes the frozen product canister for manual transfer for drying in a traditional lab freeze dryer

The UI has been designed to be intuitive not only in terms of the hardware interaction for first time users, but also thanks to its multi-touch HMI and integrated LED messaging board on the front panel.

Integrated instrumentation provides the user with real-time process data, including the mass of the liquid product and temperature of the temperature controlled liquid product reservoir platform.

Built in connectivity features allow for remote support from IMA's global service team.



LYnfinity

THANK YOU

dennis.posheluk@ima.it