

# **Advancements in Isolator** Technology

Julien Cabon **Business Development Manager** 

# Agenda

- Gloveless/Robotic Isolators
- RFID Technology
- Biological Indicator Alternatives
- Polypropylene construction
- Particle monitoring
- Ingress/Penetration testing





# Gloveless/Robotic Isolator Systems



## **Challenges with gloved or open isolators**

- Biggest risk to an aseptic process is... the human!
- Standard isolators carry a risk of gloves rupturing, contaminating the process
- Routine glove testing is
  required which takes time





# **Separation Assurance Continuum**

Increasing separation increases assurance of sterility in process





LIFE SCIENCES

Theoretical robot in an autoclave



# **Increased efficiency**

- Can perform a repetitive process quicker than a human, increasing capacity
- Can reduce human mistakes in critical processes e.g. fill/finish
- But are still prone to breakdowns and process mistakes
- Robotic isolators eliminate the human operator from the process, minimising risk
- Contamination risk still remains from interventions i.e. loading of materials
- Robotic loading of isolator is the next step to further remove risk from process



robot inside a stainless steel isolator



### LIFE SCIENCES

# Example – Aseptic Technologies Crystal L1

# **Industry Uptake**

- Five contract manufacturing organisations (CMOs) have publicly announced their use of Vanrx's SA25 Aseptic Filling Workcell.
  - **AB Biotechnologies** •
  - Fujifilm Diosynth •
  - Patheon (Thermofisher) •
  - **Singota Solutions**
  - Wuxi Biologics •





# 

# **RFID Technology**



### **Radio Frequency Identification Technology**

 $\diamond$  Digital data encoded in tags on H<sub>2</sub>O<sub>2</sub> bottles captured by a reader on isolator



Benefits include:

- Prevents bio-decontamination cycle from starting without a sufficient quantity of  $H_2O_2$  solution
- Prevents expired peroxide from being used
- Traceability by printing batch number on cycle report
- Prevents other potentially dangerous chemicals from being accidentally vaporised

Ensures only approved peroxide is used...



### LIFE SCIENCES

\*\*\*\*\*\* \*\*\*\*\*\* L-4 Opti \*\*\*\*\*\* 11/12/2020 09:52:02 PLC VERSION: V1.08 HMI VERSION: V2.02 SERIAL NUMBER: 201601L44172 CYCLE COUNT: 1 CALIBRATION DUE DATE: 09/12/2021 EOUIPMENT REFERENCE: ABF0950110111091 H202 LOT NUMBER: AAA0000000000000 CYCLE NAME: Training test cycle 3. PRESSURE TEST ENABLED 4. TEST PRESSURE 100 Pa TEST TIME 180 s TEST DECAY 50 Pa AIRFLOW SETPOINT 20 m<sup>3</sup>/h .DELIVERY TEMPERATURE 60°C

10.PRESSURE CONTROL ENABLEI PRESSURE LOW ALARM 30 Pa PRESSURE HIGH ALARM 50 Pa CONDITIONING TIME 300 s GASSING INJECTION RATE 2.0 g/m .GASSING TIME 180 s 19.DWELL INJECTION RATE 1.5 g/m 20.DWELL TIME 180 s 21.H202 ALERT LEVEL 24.AERATION VFC DELAY 300 s 25.AERATION PARAMETRIC ENABLED 26.AERATION TIME 600 s

### Isolator bottle module

USER: Engineer

SIGNATURI

\*\*\*\* 09:52:14 PRESSURISING 3 Pa \*\*\*\*

10:05:53 TEST START 100 Pa 

10:14:30 GASSING 41Pa 20.9m<sup>3</sup>/hr 57.5°c 49%

0.0ppm 

10:17:30 GASSING DWELL 39Pa 20.0m<sup>3</sup>/hr 61.1°c 85% mqq0.0

Example report

### **Regulatory requirements - BPR**

EU regulation 528/2012 Biocidal Products Regulation (BPR) came into force in Europe in September 2013 and applies to all biocidal products used in Europe

- The regulation concerns both the sale and use of biocidal products in Europe and involves: The analysis of a product's performance (efficacy)
  - Toxicity
  - Environmental fate and risk during use
- Registrations include not only the biocide, but the method in which it is introduced to the enclosure
- Therefore RFID technology ensures that only the registered biocide is used in combination with the registered equipment/bio-decontamination method





# **Biological Indicators** Alternatives

# **Chemical Indicators**

Why do we use chemical indicators?

Their use is recommended in regulations

### **FDA Guidance for Industry states:**

"An appropriate, quantified Biological Indicator (BI) challenge should be placed on various materials"

"Chemical indicators may also be useful as a qualitative tool to show that the decontaminating agent reached a given location. "

Calibrated Chemical Indicators are a new technology which are a quantitative tool...



### SCIENCES

# **Calibrated Chemical Indicators**

 Consist of a card indicator containing a chemical strip impregnated with ink

Ink reacts to hydrogen peroxide vapour and changes colour to indicate the efficacy of the biodecontamination cycle.

 Provides instant indication of the biodecontamination process efficacy

Results are calibrated against biological indicators





# How do they work?

### Reactive ink

- The reactive ink portion changes colour
- It will all change to the same colour at the end
- Reference colours
  - The reference colours do not change colour
  - ✓ The number "6" indicates log 6
  - The number ">6" indicates beyond log 6 was reached
  - The numbers "2" and "4" indicate log 6 is not reached







# **Chemical indicator reader**

- A handheld device which enables customers to quickly and easily obtain a real-time, quantifiable log reduction indication from a biodecontamination cycle
- Scans reactive ink and uses software to calculate equivalent BI log
- Provides a PDF report of the results including: •
  - Equivalent log reduction
  - Date
  - User login
  - A copy of the instrument calibration certificate





**CI** reader

Example report

**Bioquell** CI Report

User name: Doe, John Device Calibration due: 26/04/2023

**Bioquell Room Cl** Lot number: 12345L1234

Name	Date	Time	dE	Log reduction	Result
Sample 1	08/04/2019	15:23	10.18	7.27	PASS
Sample 2	08/04/2019	15:23	8.98	7.46	PASS
Sample 3	08/04/2019	15:24	9.43	7.39	PASS
Sample 4	08/04/2019	15:24	9.47	7.38	PASS
Sample 5	08/04/2019	15:24	7.60	7.69	PASS
Sample 6	08/04/2019	15:24	8.21	7.59	PASS
Sample 7	08/04/2019	15:25	8.82	7.49	PASS
- · · ·	001010010		0.00		****

### LIFE SCIENCES

### Report Date: 03/02/2022 Device S/N: 2010006614

### Expiry date: 25/11/2022

# **Enzymatic Indicators**

- Based on the inactivation of a tAK (thermostable Adenylate Kinase) enzyme
- Enzyme catalyses a reaction that produces light and the amount of light produced is directly comparable to the amount of remaining tAK activity on the indicator.
- Light output is measured using a luminometer and provides a quantifiable numerical value
- Results are outputted to a computer which correlates light units to the log reduction of a BI



# **Comparison between indicators**

Category	BI	Calibrated CI		
Mode of activity	Inactivation of biological spore	Breakdown of chemical dye		
Time to result	7 days	Instant Read		
Read method	By eye after incubation in TSB at 57°C	By eye or by light reader		
Quantifiable result	No	Yes (vs. BI)		
Shelf Life	12 months	12 months		
Use supported by Regulations?	Yes	No		
Cost per indicator	\$10	\$2		



### LIFE SCIENCES

EI

### Inactivation of enzyme

# Minutes (after 30 min test prep)

React with Luciferin and read with luminometer

Yes

6 months

No

\$25 (plus reader @ \$20k)



**Polypropylene construction** 

# Background

- Traditionally, all isolator systems were made of stainless steel
- Developments in technology have led some companies to develop isolators from a unique polypropylene construction
- This has several advantages over • stainless steel...





# **Stainless Steel vs Polypropylene comparison**

### **Stainless Steel**

### Polypropylene

### **Manufacturing process**

Fabricated isolators will all be slightly different even if manufactured to the Polypropylene isolators are manufactured from the same mould, therefore, same design as they are 'handmade' – cut, bent, welded, polished etc, which the final product is consistent and correct every time. However, this means can result in

inconsistent fit and variation between each unit

# the dimensions will be standard, and customization is more complex

### **Manufacturing time**

Stainless Steel isolators require significant lead times based on the production/fabrication process. Fabrication alone can take up to 6 weeks

Polypropylene isolators can be manufactured in as little as one day allowing delivery in as little as 6 weeks from the time of order through validation





# **Stainless Steel vs Polypropylene comparison**

### **Stainless Steel** Polypropylene Weight Stainless steel is extremely heavy putting significant load on floor Polypropylene has a density 8x less than stainless steel, so is optimal for stainless steel, so is optimal for minimal floor loads loads **Ergonomics** Standard box shapes mean corners can be sharp and can allow for trapped The moulding process allows smooth corners, which makes cleaning easier and can allow for trapped dirt. which makes cleaning easier (GMP compliant) Leak tightness Each weld, joint and gasket can lead to a potential risk of leaks A single-piece moulded chamber shell means there are fewer penetrations, so risk of leaks is drastically reduced





# ECELAB® 6 bloquel An Ecolab Solution

# **Particle monitoring**



# **Particle Monitoring – Regulatory changes**

Updated Annex 1 section 5.9 states sampling head must be "as close as possible to critical location"

5.9 Particle counters, including sampling tubing, should be qualified. The manufacturer's recommended specifications should be considered for tube diameter and bend radii. Tube length should typically be no longer than 1m unless justified and the number of bends should be minimized. Portable particle counters with a short length of sample tubing should be used for classification purposes. Isokinetic sampling heads should be used in unidirectional airflow systems. They should be oriented appropriately and positioned as close as possible to the critical location to ensure that samples are representative.

### However, Annex 1 section 4.19 also states isolators should ensure "first air protection in the critical zone"

4.19 The design of the technology and processes used should ensure appropriate conditions are maintained in the critical zone to protect the exposed product during operations.

- Isolators: 1.
  - a. The design of open isolators should ensure grade A conditions with first air protection in the critical zone and unidirectional airflow that sweeps over and away from exposed products during processing.



# **Centrally Located Isokinetic Cone**

- To comply with 4.19 and 5.9 of Annex 1, some isolator manufacturers have developed a centrally located iso cone position
- Consists of a sampling pipe with bends
- Position can be adjusted during install/OQ to allow operators to define sampling location, but must be fixed during operation.
- Location should be a compromise between being as close as possible to 'critical zone' without compromising 'first air'





# 

# Ingress/Penetration testing

# Background

EU GMP Annex 1 section 4.22 requires customers to demonstrate that bio-decontamination agents do not impact the medicinal product i.e. ingress into packaging

> 4.22 Decontamination methods (cleaning and bio-decontamination, and where applicable inactivation for biological materials) should be appropriately defined and controlled. The cleaning process prior to the bio-decontamination step is essential; any residues that remain may inhibit the effectiveness of the decontamination process. Evidence should also be available to demonstrate that the cleaning and biodecontamination agents used do not have adverse impact on the product produced within the RABS or isolator.

 $A_2O_2$  vapour is the most common bio-decontamination agent used in isolators

 $\diamond$ Industry requires an H<sub>2</sub>O<sub>2</sub> vapour penetration testing service...



# **Hydrogen Peroxide Vapor Ingress Testing**

Determines concentration of  $H_2O_2$  (down to ppb) which may permeate into critical samples i.e. vials and bags, which are decontaminated with  $H_2O_2$  vapor

Sample is subjected to  $H_2O_2$  vapor biodecontamination and a horseradish peroxidase assay is performed to determine the levels of  $H_2O_2$  present inside sample





# Hydrogen Peroxide Vapor Ingress Testing

	Test required		H <sub>2</sub> O <sub>2</sub> Amplex r (developed) or	
Report is provided to summarize what parameters were used in the bio-decontamination cycle, and what concentration of $H_2O_2$ (if any) was measured inside the sample tested	Sample Extract Reading (ppm)	5ml Test Gass Un-g Test Gass Un-g	5ml Ampoule: Test 1, 25-June-2013 Gassed; 0.043, 0.039 Un-gassed; 0.010, 0. Test 2, 27-June-2013 Gassed; 0.031, 0.031 Un-gassed; 0.007, 0. Test 3, 1-July-2013; Gassed; 0.041, 0.047 Un-gassed; 0.005, 0. Notes: Overall average of ga All un-gassed sample (0.015ppm) Assay detection limits	
Customers then decide whether the exposure level of $H_2O_2$ is acceptable		Test Gass Un-g Note Over All u (0.0		





### LIFE SCIENCES

red assay test on sample following HPV gassing cycle n sample within Walk-in-Chamber (WiC)

9, 0.040 (ppm) (average of 0.041ppm) 010, 0.009 (ppm)

1, 0.031 (ppm) (average of 0.031ppm) .007, 0.007 (ppm)

7, 0.043 (ppm) (average of 0.044ppm) 004, 0.005 (ppm)

ssed samples = 0.039ppm s have a reading below the lower detectable limit

s = 0.015 - 0.189 ppm.

# **Example report**

# Thank you and any questions?



